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Instructions for Use – Please read carefully

Abutments
Adapters, screws and gingiva formers

NOTE: For the stabilisation of removable dentures at least 6 implants per jaw are necessary. If less implants are used, abutments may show irreparable damage through wear.

I. Indications for Use
Abutments and prosthetic screws and other components are intended for use with dental implants as a support for single or multiple tooth prostheses in the maxilla or mandible of a partially or fully edentulous patient.

II. Contraindications
The abutments of each series can only be combined with the corresponding compatible implant system. No inappropriate in connection geometry abutments should be used. Any post-processing at the connection geometry to the implant may result in fitting inaccuracies prohibiting further use or to damages. Contraindicated in adolescents and children with a not fully developed jaw.

III. How Provided
Abutments are provided with the abutment screw for attachment to the dental implant (where relevant).

IV. Composition
Allfit® Abutments and abutment screw: Ti-6Al-4V ELI, AF1M F 136.

V. Warning
Abutments are intended for single use. Side effects: Allergies or sensitibilities in connection with the alloy are rare but cannot be ruled out. Interactions: Various alloy types in the same mouth cavity can lead to galvanic reactions in the event of occlusal or proximal contact.

VI. Manufacturer's Safety Statement
Abutments have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the abutment in the implant's environment is unknown. Scanning a patient who has this device may result in patient injury.

VII. Torque
Only the implant manufacturer's recommended torque is to be used.

Ncm	Implant Diameter
15	≥ 3,0 mm
20	3,3 mm
25	3,5 - 4,8 mm
30	5,0, 5,5 mm

VIII. Application
Caution: The following descriptions are not sufficient for the immediate use of the abutment and other Suprastructures. Dental skills and prior instruction of how to use the abutments are required.
For the attachment of overdenture or partial dentures on endosseous dental implants: The abutments are directly screwed into the implant and have a defined mounting various matrix systems. These are attached to the abutments using the matrix systems that are inserted in the overdentures or partial dentures. Every type of reworking of the connection geometry to the implant leads to inaccurate fits which rule out re-use. In addition the facing defined connection points to the matrix systems may not be ground or changed as they are imperative to guarantee a sufficient amount of stability and a good fit. For the stabilisation of removable dentures at least six implants per jaw are necessary. If ten implants are used, abutments may show a maximum of durability.

IX. Cleaning and Sterilization
The abutments are delivered in a non-sterile state in the appropriate packaging. The product must be cleaned and sterilized according to following protocol:

Cleaning:
Automatic cleaning procedure based on Vario TD program: Cold water rinse 1 min. Cleaning at 55° C/131° F (± 2° C/35,6° F) for 5 min. With alkaline cleaning agent, neutralization with cold water for 2 min. Cold water rinse 1 min.

Sterilization:
Seal in a suitable autoclaveable pouch. Steam sterilization with 3x fractionated pre-vacuum: exposure at 134-137° C (270° F) for 4 minutes, drying time: min. 15 minutes.

X. Precautions
The product must be inspected prior to usage. The packing must be sealed without any visible damage. The following descriptions are not sufficient for the immediate use of the product. Dental skills and prior instruction of how to use the product are at any rate required.

Legend

	Read instructions
	Expiration date
	Only use once
	non sterile
	This product is only for sale to dentists, oral surgeons and maxillo-facial surgeons or on their behalf.
	Charge number
	Store in dry place
	Manufacturer
	Catalogue number
	-5° - 55°
	Authorized representative in the European Community

XI. Notice
The information conveyed orally, in writing or in practical seminars is based on tests and experience and can therefore only be considered to be standard values. Our products are subject to constant further development. In this connection we reserve the right to make product changes in relation to the design and composition.

Keep abutments out of reach of children. To be used by dentists or surgeons only.

996-ABUTMENTS-03
Revision issued 15-03-2022

SIMPLADENT GmbH
Dorfplatz 11
CH-8737 Gommiswald
SWITZERLAND
Tel. +41 (0)55 293 23 70
Fax +41 (0)55 293 23 00
simpladent@implant.com
www.simpladent-implant.com

Distributed by
Simpladent Implant Solutions Pvt Ltd
Unit no 417, 4th Floor, KM Trade Tower,
H-3, Sector 14, kaushambi, Ghaziabad
201010 Uttar Pradesh
India

EC	REP	Dr. Ihde Dental GmbH Erfurter Str. 19 D - 85386 Eching / München Tel +49 (0)89 319 761 0 info@ihde-dental.de
-----------	------------	--